

Applied Orthogonal Experiment Design for the Optimum Extraction Conditions of High Concentration Selenium from Maifanite

Hailin Liao^{1,2}, Lingfeng Jiang^{1,2*}, Chunjie Wei^{1,2}, Hua Huang^{1,2}, Jing Pan^{1,2}, Chuntao Luo^{1,2}

¹Guangxi Scientific Research Centre of Traditional Chinese Medicine, Guangxi University of Chinese Medicine, Nanning, China

²Key Laboratory of Generic Technology Research and Development of Traditional Chinese Drugs Pharmaceutics in Guangxi Colleges and Universities, Guangxi University of Chinese Medicine, Nanning, China

Email: *jianglingfeng@126.com

How to cite this paper: Liao, H.L., Jiang, L.F., Wei, C.J., Huang, H., Pan, J. and Luo, C.T. (2017) Applied Orthogonal Experiment Design for the Optimum Extraction Conditions of High Concentration Selenium from Maifanite. *Detection*, 5, 1-4.
<https://doi.org/10.4236/detection.2017.51001>

Received: January 11, 2017

Accepted: January 28, 2017

Published: January 31, 2017

Copyright © 2017 by authors and Scientific Research Publishing Inc.

This work is licensed under the Creative Commons Attribution International License (CC BY 4.0).

<http://creativecommons.org/licenses/by/4.0/>



Open Access

Abstract

In order to improve the concentration of selenium in the maifanite concentrate, the optimum experimental conditions were obtained by orthogonal experiment as follows: size of maifanite was 10 mesh, concentration of maifanite was 15 (expressed as the quality of maifanite divided by that of water), soaking time was 60 min at room temperature, heating temperature was 70 °C, heating time was 20 min, pH value was 6. Under this condition, the concentrations of each element (ppm) determined by ICP-AES were: Fe 0.099, Cu 0.035, Mn 0.051, Zn 0.019, Se 0.028, Al 0.000. a ven in this document.

Keywords

Selenium, Maifanite, ICP-AES, Style, Styling

1. Introduction

Metal dyshomeostasis was one of pathological characterizations of Alzheimer's disease (AD) [1] [2] [3] [4]. Maifanite, a mineral medicine and was used in many fields for a long time in East Asia countries such as China and Japan [5] [6]. It might be the ideal metal homeostasis reagent. In our previous study, we found that maifanite could balance the dyshomeostasis of trace elements in A β -induced AD rats, and the trace element Se was very important, and its level in high dose maifanite group was higher than that in normal control group [7]. There has been heightened interest in the role of the trace element Se in health and neurologic disorders including AD [1].

Many mineral elements were produced after the primary maifanite mineral was immersed into water [6]. The major elements in the aqueous extract of mai-

fanite were K, Na, Ca, Mg, Sr, Zn, Cu, Fe, Mn and Se etc. [4] [5] [6]. In order to improve the concentration of selenium in the maifanite concentrate, we designed the orthogonal experiments from six factors: size of maifanite, concentration of maifanite, soaking time at room temperature, heating temperature, heating time and pH value. The optimum dissolution conditions of maifanite were studied by using Se content as the index of dissolution performance.

2. Experiment

The maifanite was purchased from the Guifeng Maifanite Company in Guangxi Guiping City, China. Deionized water (18 M Ω cm, Milli-Q gradient system, Millipore) was used throughout the experiments.

By the inductively coupled plasma-atomic emission spectroscopy (ICP-AES; IRIS Intrepid II XSP, USA Thermo Elemental), the elements Fe, Cu, Zn, Mn, Se, and Al were analyzed in the optimum dissolution conditions of maifanite. The measurement conditions were 1200 W from the radio-frequency generator, a plasma argon flow rate of 15 L/min, a cooling gas flow of 14 L/min, a carrier gas flow of 1.0 L/min, a 20- μ m entrance slit, a 30- μ m exit slit, a height of observation of 15mm, and an integration time lapse of 5s.

Software of orthogonal designing assistant II V3.1 was used for the evaluation of the statistical experimental design. A five-level and six-factor orthogonal experiment (**Table 1**) was designed according to the L25-5-6 table as follow.

3. Results

According to the orthogonal experiment design (L25-5-6 table), twenty-five experiments should be accomplished. As shown in **Table 2**, the results of Se content were determined by ICP-AES in water extract of maifanite. Regarded Se content as the index of dissolution performance, the optimum dissolution conditions of maifanite were obtained as follow: size of maifanite was 10 mesh, concentration of maifanite was 15 (expressed as the quality of maifanite divided by that of water), soaking time was 60 min at room temperature, heating temperature was 70°C, heating time was 20 min, pH value was 6. Under this condition, the concentrations of each element (ppm) determined by ICP-AES were: Fe

Table 1. A five-level and six-factor orthogonal experiment was designed according to the L25-5-6 table by orthogonal designing assistant II V3.1.

No.	1	2	3	4	5	6
Factor name	Size/mesh	Concentration ¹⁾	soaking time/min ²⁾	heating temperature/°C	heating time/min	pH value
level 1	10	1	10	50	20	5.00
level 2	50	5	30	60	40	6.00
level 3	100	10	60	70	60	7.00
level 4	150	15	90	80	80	8.00
level 5	200	25	120	90	120	9.00

¹⁾expressed as the quality of maifanite divided by that of water. ²⁾soaking time at room temperature.

Table 2. The results of Se content were determined by ICP-AES in water extract of maifanite.

column	1	2	3	4	5	6	
Factor name	Size/mesh	Concentration ¹⁾	soaking time/min ²⁾	heating temperature/°C	heating time/min	pH value	Results/ppm
experiment 1	1	1	1	1	1	1	0.035
experiment 2	1	2	2	2	2	2	0.022
experiment 3	1	3	3	3	3	3	0.038
experiment 4	1	4	4	4	4	4	0.036
experiment 5	1	5	5	5	5	5	0.014
experiment 6	2	1	2	3	4	5	0.031
experiment 7	2	2	3	4	5	1	0.023
experiment 8	2	3	4	5	1	2	0.026
experiment 9	2	4	5	1	2	3	0.027
experiment 10	2	5	1	2	3	4	0.025
experiment 11	3	1	3	5	2	4	0.022
experiment 12	3	2	4	1	3	5	0.019
experiment 13	3	3	5	2	4	1	0.021
experiment 14	3	4	1	3	5	2	0.019
experiment 15	3	5	2	4	1	3	0.015
experiment 16	4	1	4	2	5	3	0.027
experiment 17	4	2	5	3	1	4	0.024
experiment 18	4	3	1	4	2	5	0.022
experiment 19	4	4	2	5	3	1	0.029
experiment 20	4	5	3	1	4	3	0.027
experiment 21	5	1	5	4	3	2	0.023
experiment 22	5	2	1	5	4	3	0.026
experiment 23	5	3	2	1	5	4	0.017
experiment 24	5	4	3	2	1	5	0.038
experiment 25	5	5	4	3	2	1	0.029
mean value 1	0.029	0.028	0.025	0.025	0.028	0.027	
mean value 2	0.026	0.023	0.023	0.027	0.024	0.018	
mean value 3	0.019	0.025	0.030	0.028	0.027	0.032	
mean value 4	0.026	0.030	0.027	0.024	0.027	0.025	
mean value 5	0.027	0.022	0.022	0.023	0.020	0.025	
range	0.010	0.008	0.008	0.005	0.008	0.014	

¹⁾expressed as the quality of maifanite divided by that of water. ²⁾soaking time at room temperature.

0.099, Cu 0.035, Mn 0.051, Zn 0.019, Se 0.028, Al 0.000.

4. Conclusion

According to the orthogonal experiment design, regarded Se content as the index of dissolution performance, the optimum dissolution conditions of maifanite were obtained as follows: size of maifanite was 10 mesh, concentration of maifanite was 15 (expressed as the quality of maifanite divided by that of water), soaking time was 60 min at room temperature, heating temperature was 70 °C,

heating time was 20 min, pH value was 6. Under this condition, the concentrations of each element (ppm) determined by ICP-AES were: Fe 0.099, Cu 0.035, Mn 0.051, Zn 0.019, Se 0.028, Al 0.000 ure/K”.

Acknowledgements

This work was supported by the Sub Project of the Key Laboratory of Generic Technology Research and development of Traditional Chinese Drugs Pharmaceuticals in Guangxi colleges and Universities (No. ZJGX201402005), which affiliated to Promotion Project of Science and Technology Innovation Ability in Guangxi Colleges and Universities (approved by the Education Department of Guangxi Government, China, No. [2014] 14).

References

- [1] Loef, M., Schrauzer, G.N. and Walach, H. (2011) Selenium and Alzheimer’s Disease: A Systematic Review. *Journal of Alzheimers Disease*, **26**, 81-104.
- [2] Ozcelik, D., Uzun, H. and Naziroglu, M. (2012) N-Acetylcysteine Attenuates Copper Overload-Induced Oxidative Injury in Brain of Rat. *Biological Trace Element Research*, **147**, 292-298. <https://doi.org/10.1007/s12011-012-9320-1>
- [3] Takahashi, S., Takahashi, I., Sato, H., Kubota, Y., Yoshida, S. and Muramatsu, Y. (2001) Age-Related Changes in the Concentrations of Major and Trace Elements in the Brain of Rats and Mice. *Biological Trace Element Research*, **80**, 145-158. <https://doi.org/10.1385/BTER:80:2:145>
- [4] Koç, E.R., İlhan, A., Zübeyde, A., Acar, B., Gürler, M., Altuntaş, A., Karapirli, M. and Bodur, A.S. (2015) A Comparison of Hair and Serum Trace Elements in Patients with Alzheimer Disease and Healthy Participants. *Turkish Journal of Medical Sciences*, **45**, 1034-1039. <https://doi.org/10.3906/sag-1407-67>
- [5] Juan, L., Zhang, P.Y., Gao, Y., Song, X.G. and Dong, J.H. (2008) Overview of Maifanshi: Its Physi-Chemical Properties and Nutritious Function in Drinking Water. *Environmental Science & Technology*, **31**, 63-66.
- [6] Zhang, B.G. (2005) The Progress of Pharmaceutical Research on the Maifanite. *Chinese Traditional Patent Medicine*, **27**, 1205-1208.
- [7] Jiang, L.F., Liao, H.L., Huang, H.M., Zhou, L.X., Li, L., Cheng, S.X. and Du, C.Z. (2013) Potential Prevention and Treatment of Maifanite for Alzheimer’s Disease Based on Behavior Test, Oxidative Stress Assay and Trace Element Analysis in Hippocampus of Abeta (25-35)-Induced AD Rats. *Biological Trace Element Research*, **152**, 50-56. <https://doi.org/10.1007/s12011-012-9590-7>

Submit or recommend next manuscript to SCIRP and we will provide best service for you:

Accepting pre-submission inquiries through Email, Facebook, LinkedIn, Twitter, etc.

A wide selection of journals (inclusive of 9 subjects, more than 200 journals)

Providing 24-hour high-quality service

User-friendly online submission system

Fair and swift peer-review system

Efficient typesetting and proofreading procedure

Display of the result of downloads and visits, as well as the number of cited articles

Maximum dissemination of your research work

Submit your manuscript at: <http://papersubmission.scirp.org/>

Or contact detection@scirp.org